RESEARCH ARTICLE

Open Access

Fully-automated production of [⁶⁸Ga]Ga-FAPI-46 for clinical application



Sarah Spreckelmeyer * Matthias Balzer, Simon Poetzsch and Winfried Brenner

* Correspondence: sarah. spreckelmeyer@charite.de Charité - Universitätsmedizin Berlin, corporate member of Freie Universität Berlin, Humboldt-Universität zu Berlin, and Berlin Institute of Health, Department of Nuclear Medicine, Augustenburger Platz 1, 13353 Berlin, Germany

Abstract

Background: [⁶⁸Ga]Ga-FAPI-46 is a promising radiopharmaceutical for in vivo detection of the fibroblast activation protein by positron emission tomography. Until now, the synthesis of [⁶⁸Ga]Ga-FAPI-46 has been only performed manually. Our aim was to evaluate the automated synthesis of this radiopharmaceutical on two different commercially available synthesis modules in order to make the tracer readily available for clinical application.

Results: The synthesis of [⁶⁸Ga]Ga-FAPI-46 with different amounts of precursor (10–50 µg) on the Modular Lab PharmTracer (MLPT) and Modular Lab eazy (ML eazy) from Eckert & Ziegler with a customized synthesis template and a customized single-use cassette yielded best results with 50 µg FAPI-46 for clinical multi-dose application. All relevant quality control parameters tested (e.g. sterility, stability and radiochemical purity) were in accordance with the European Pharmacopoeia.

Conclusions: [⁶⁸Ga]Ga-FAPI-46 was successfully synthesized fully-automated on the synthesis modules Modular Lab PharmTracer and ML eazy and is, thus, available for multi-dose application in clinical settings.

Keywords: Automated production, Imaging, FAP, FAPI-46, Fibroblast

Background

The fibroblast activation protein (FAP) is an endopeptidase and is located in the cell membrane. FAP is not expressed in normal, healthy tissue, but during wound repair and in the microenvironment of numerous cancer types, including colorectal and pancreatic cancer (Brennen et al. 2012). It has been shown, that FAP is selectively targeted by fibroblast activation protein inhibitors (FAPIs) with a *N*-(4-quinolinoyl)glycyl-(2-cyanopyrrolidine) scaffold (Jansen et al. 2013). In 2018, Lindner et al. developed a series of quinolone-based molecules (e.g. FAPI-04, Fig. 1) for FAP-based imaging and therapy with gallium-68 and lutetium-177, respectively (Lindner et al. 2018). One year later, the same group optimized the structure activity relationships (SARs) and succeeded to improve FAP binding and pharmacokinetics yielding high-contrast images. The most promising molecule from that study is FAPI-46 (Fig. 1) (Loktev et al. 2019). Recently, Giesel et al. reported of FAPI-74, which main advantage is that it can be label



© The Author(s). 2020 **Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit http://creativecommons.org/licenses/by/4.0/.

with gallium-68 and fluor-18, resulting in [⁶⁸Ga]Ga-FAPI-74 and [¹⁸F]F-FAPI-74, respectively (Giesel et al. 2020).

To our knowledge, the synthesis of [⁶⁸Ga]Ga-FAPI-46 has not yet been described on an automated synthesis module. Here, we would like to introduce the synthesis conditions of [⁶⁸Ga]Ga-FAPI-46 on two different synthesis modules, namely Modular Lab PharmTracer (MLPT) and Modular Lab eazy (ML eazy) from Eckert & Ziegler Eurotope GmbH – in order to introduce a fully automated [⁶⁸Ga]Ga-FAPI-46 synthesis for clinical application. For that purpose, different loads of FAPI-46 were evaluated based on quality control parameters and the radioactivity distribution on the different cassette parts after synthesis. For the synthesis of ⁶⁸Ga-tracers, different modules are available on the market, which have been reviewed by Boschi et al. (Boschi et al. 2013) In addition, generator post-processing steps like fractionation, anionic-exchange and cation-exchange are described (Boschi et al. 2013). The now proposed fully automated pre-purification and concentration of gallium-68 on a strong cation exchange (SCX) cartridge has the advantage that gallium-68 becomes concentrated while germanium-68 as well as non-radioactive impurities (e.g. Zn²⁺) are trapped on the SCX cartridge (Velikyan 2015).

Results

Labeling results of different FAPI-46 loads on MLPT

The fully automated production of $[^{68}$ Ga]Ga-FAPI-46 was conducted on the commercial labeling synthesis modules MLPT and ML eazy. Five different precursor amounts of FAPI-46 were investigated for radiolabeling with gallium-68 (10–50 µg, n=3) on MLPT. Immediately after synthesis, the different cassette parts were analyzed for radioactivity in order to get an overview of the radioactivity distribution on the cassette and reconcile these findings with the effects of the labeling conditions that were used. As shown in Table 1, 30 µg, 40 µg and 50 µg of FAPI-46 resulted in a similar radioactivity distribution in the product vial $-95.1\pm0.4\%$, $95.3\pm1.2\%$ and $96.0\pm0.6\%$, respectively. On the CM cartridge, the radioactivity decreases from $2.9\pm0.3\%$ to $2.1\pm0.8\%$. By lowering the amount of FAPI-46 to $10\,\mu$ g, the percentage found in the product decreased significantly to $70.2\pm6.1\%$, and the radioactivity on the CM cartridge increased to $26.2\pm6.1\%$. The radioactivity found in the reaction vial, in the waste fraction and on the SCX cartridge however did not change significantly by changing the amount of the starting material.

Labeling results of different FAPI-46 loads on ML eazy

On the ML eazy module, we performed experiments with $10-50 \,\mu g$ of FAPI-46 (n=3) at 98 °C reaction temperature (Table 2). As a result, $20-50 \,\mu g$ FAPI-46 load show the

Table 1 Activity measured on different parts of the cassette immediately after synthesis on MLPT (n = 3)

load [μg]	product vial [%]	CM cartridge [%]	waste [%]	SCX cartridge [%]	reaction vial [%]	RCY [%]
10	70.2 ± 6.1	26.2 ± 6.1	1.5 ± 1.0	0.6 ± 0.1	1.5 ± 0.3	68.9 ± 3.8
20	91.5 ± 1.9	4.8 ± 0.8	2.1 ± 1.5	0.6 ± 0.3	1.0 ± 0.3	91.8 ± 3.9
30	95.1 ± 0.4	2.9 ± 0.3	0.4 ± 0.2	0.4 ± 0.1	1.1 ± 0.4	94.2 ± 0.8
40	95.3 ± 1.2	2.7 ± 0.9	0.6 ± 0.1	0.3 ± 0.0	1.1 ± 0.3	95.1 ± 2.1
50	96.0 ± 0.6	2.1 ± 0.8	0.5 ± 0.1	0.5 ± 0.2	0.9 ± 0.2	95.2 ± 1.4

highest radioactivity distribution in the product vial. Decreasing the FAPI-46 load to 10 µg leads to a significant decrease in radioactivity measured in the product vial.

Comparing the distribution pattern of the MLPT with ML eazy (Fig. 2), the results are similar with $20-50\,\mu g$ FAPI-46 load. Whereas with $10\,\mu g$, a significant difference can be observed ($70.2\% \pm 6.1\%$ on MLPT vs. $36.1\% \pm 15.5\%$ on ML eazy).

Quality control of [68Ga]Ga-FAPI-46

The radiochemical purity was evaluated by radio-HPLC and radio-iTLC.

With *radio*-HPLC, free gallium-68 was detected at t_R = 3.0 min, whereas gallium-68 bound to FAPI-46 was detected at t_R = 6.9 min with a purity of 99.7%. Radioactive impurities were detected at t_R = 5.0 min with approx. 0.1% and t_R = 6.6 min with approx. 0.2% of the total radioactivity. (Fig. 3a).

The cold standard [nat Ga]Ga-FAPI-46 shows a similar retention time of $t_R = 6.5$ min as seen in Fig. 3b.

With *radio*-iTLC, no free gallium-68 (Fig. 4a) or 68 Ga-colloide (Fig. 4b) could be detected at R_f = 0.8 or R_f = 0.2, respectively. The product was detected at R_f = 0.2 (Fig. 4a) and R_f = 0.8 (Fig. 4b).

In addition, the product solution was tested for endotoxins. For this purpose, the solution was diluted with endotoxin-free-water in a ratio 1:100. For all samples, the endotoxin concentration was below 5.0 IE/mL. This is in accordance with the European Pharmacopoeia (9.0/0125), where a threshold of 175 IE/mL is given. With regard to sterility, all products were sterile.

The stability of [⁶⁸Ga]Ga-FAPI-46 in aqueous solution at room temperature was tested up to 3 h by *radio*-HPLC. As seen in Fig. 5, [⁶⁸Ga]Ga-FAPI-46 is stable under those conditions. No additional radioactive by-products or free gallium-68 could be detected during this time period.

Table 2 Activity measured on different parts of the cassette immediately after synthesis on ML eazy at $98 \,^{\circ}\text{C}$ (n = 3)

load [µg]	product vial [%]	CM cartridge [%]	waste [%]	SCX cartridge [%]	reaction vial [%]	
10	36.1 ± 15.5	59.7 ± 15.6	0.0 ± 0.0	0.4 ± 0.0	3.7 ± 0.2	
20	75.3 ± 18.2	19.3 ± 14.2	0.2 ± 0.2	0.4 ± 0.1	4.8 ± 4.2	
30	88.7 ± 10.1	9.2 ± 9.7	0.0 ± 0.0	0.5 ± 0.1	1.6 ± 0.5	
40	85.5 ± 11.9	11.0 ± 9.6	0.1 ± 0.1	0.5 ± 0.2	2.9 ± 3.0	
50	89.7 ± 6.7	6.1 ± 5.5	0.0 ± 0.0	0.5 ± 0.1	3.7 ± 1.4	






Fig. 2 Radioactivity distribution on cassette parts after synthesis on MLPT (95 °C) and ML eazy (98 °C) (dark: MLPT; light: ML eazy) with increasing peptide load from left to right (10 μ g–50 μ g); *= p-value < 0.05

The product specifications of $50 \,\mu g$ FAPI-46 are summarized in Table 3. Acceptance criteria are based on the European Pharmacopoeia (9.0/0125).


Discussion

For establishing the automated synthesis of [68 Ga]Ga-FAPI-46 on the Modular Lab PharmTracer and ML eazy from Eckert & Ziegler Eurotope GmbH, five different amounts of FAPI-46 were evaluated for the production of this tracer. The results summarized in Table 1 demonstrate that an amount of $30-50\,\mu g$ of FAPI-46 yield a similar radioactivity distribution on the cassette on MLPT, with $50\,\mu g$ yielding the lowest amount of uncomplexed gallium-68 on the CM cartridge. Thus, high quantitative radiolabeling was achieved best with $50\,\mu g$ of FAPI-46. With regard to clinical application from an economical point of view, the application of only $40\,\mu g$ or even $30\,\mu g$ FAPI-46 on MLPT looks promising and needs further evaluation by quantitative imaging data. The low standard deviations for all the productions confirm their high reproducibility.





For $10\,\mu g$ and $20\,\mu g$, an increasing amount of radioactivity on the CM cartridge and a decreasing yield of activity in the product vial was observed. Theoretically, two scenarios are possible to our knowledge for explaining the increased radioactivity on the CM cartridge. Either the product is trapped on the CM cartridge, or uncomplexed gallium-68. The product would be eluted from the CM cartridge with 0.9% NaCl solution. However, when eluting the CM cartridge manually with 5 mL of 0.9% NaCl solution after the synthesis, this resulted in no change of the amount of radioactivity measured on the CM cartridge suggesting that no product is retained on the CM cartridge. The most likely explanation therefore is that $10\,\mu g$ and $20\,\mu g$ loads are too low to complex



 $^{68}\text{Ga}^{3+}$ ions completely resulting in a trapping of unbound gallium-68 on the CM cartridge.

The minimum amount needed to complex FAPI-46 sufficiently under the described conditions was 30 μ g. In former studies in which [68 Ga]Ga-FAPI-46 was synthesized manually, 20 μ g of FAPI-46 were used with a final purification through a solid phase extraction (Oasis light HLB, Waters) and elution with ethanol (Meyer et al. 2019). The radiochemical yields were unfortunately not reported. Thus, a direct comparison to the automated process presented in this manuscript is not possible.

Comparing the radioactivity distribution of MLPT (95 °C) and ML eazy (98 °C), both modules give similar results. For the ML eazy, a higher temperature setting was needed, Since the reaction conditions were the same (e.g. buffer formulation) compared to MLPT, the reason for this finding in our opinion is related to different set-ups of the heating units of ML eazy and MLPT. The heating unit on the ML eazy develops the heat from the bottom, whereas the heating unit on the MLPT provides the heat also from the walls. Thus, we hypothesize, that a higher temperature setting is needed on the ML eazy to reach a homogenous temperature of 95 °C.

Table 3 Summary of the product specifications for 50 μ g FAPI-46, n=3

Test	Acceptance criteria	[⁶⁸ Ga]Ga-FAPI-46	
Radiochemical purity (<i>radio-HPLC</i>)	> 95%	99.7%	
Radiochemical purity (<i>radio-</i> iTLC)	> 95%	99.9%	
рН	4–8.5	4.2	
Radioactivity concentration	> 50 MBq/ mL	105–242 MBq/ mL	
Radioactivity	> 150 MBq	740–1700 MBq (depending on amount of generato used)	
Volume	2–10 mL	7 mL	
Color	Colorless	Colorless	
Molar radioactivity	1–60 MBq/ nmol	13-30 MBq/nmol	
Radionuclidic purity	> 99.9%	99.9%	
⁶⁸ Ge breakthrough	< 0.001%	0.0002%	
Endotoxins	< 19.0 IE/mL	< 5.0 IE/mL	
Stability over 3 h	> 90%	99.9%	

With regard to quality control parameters, we achieved an endotoxine-free, sterile and stable solution of [⁶⁸Ga]Ga-FAPI-46 with a radiochemical purity > 95% over a period of 3 h on both MLPT and ML eazy. Moreover, no ethanol had to be used throughout the process, which makes this synthesis favorable over the manual synthesis, because no analysis of solvent residues like ethanol by gas chromatography is necessary. Finally, all tested quality control parameters were in accordance with the European Pharmacopoeia (see Table 3).

Conclusions

[⁶⁸Ga]Ga-FAPI-46 was successfully synthesized fully-automated on the synthesis modules MLPT and ML eazy for the first time. All the tested quality parameters for radiochemical purity, pH, endotoxins and sterility were in accordance with the European Pharmacopoeia. In addition, the product solution was stable for at least 3 h after production, as shown by *radio-HPLC*. Thus, [⁶⁸Ga]Ga-FAPI-46 can be easily and reliably produced on the modules MLPT and ML eazy for routine clinical application.

Methods

Materials

FAPI-46, manufactured by ABX, and the reference standard [natGa]Ga-FAPI-46 were made available free of costs from SOFIE Biosciences. An aqueous stock solution of 1 mg/mL was prepared and kept at – 15 °C. All chemicals were of pure chemical grade, and solvents for high-pressure-liquid-chromatography (HPLC) were obtained as HPLC grade. TraceSelect water (Sigma-Aldrich) was used in all experiments. The pharmaceutical grade ⁶⁸Ge/⁶⁸Ga generator (GalliaPharm*, Eckert & Ziegler Radiopharma GmbH, Germany), Modular Lab PharmTracer (Eckert & Ziegler Eurotope GmbH, Germany) and reagent set EZ-102 (Eckert & Ziegler Eurotope GmbH, Germany) were used. The amount of detected metal impurities/ ⁶⁸Ge breakthrough as provided by the

manufacturer was less than the defined limit in the European Pharmacopeia monograph. Activity counting was performed using a borehole counter (Nuklear-Medizintechnik Dresden GmbH, Germany). HPLC was performed using the HPLC system Knauer Azura (UVD: $2.1\,L$; P6.1L) coupled with UV and radiometric (Raytest Socket 2″8103 0370) detectors. The TLC scanner MiniGita from Raytest was used. The test for endotoxins was performed using Nexgen PTS (Charles River). The CM cartridge is a cation-exchange cartridge from Sep-Pak Accell Plus light 130 mg sorbent per cartridge, $37–55\,\mu$ m particle size. The sterility tests were carried out as required by the European Pharmacopoeia by the Institut für Hygiene, Charité.

Preparation for labeling of FAPI-46 with gallium-68 on MLPT

First, the commercially available fully automated synthesis platform MLPT was equipped with a disposable single-use, customized cassette. The synthesis template and the cassette were altered as follows, as ethanol was excluded from the whole radiolabeling process. The C_{18} cartridge was exchanged by a CM cartridge. The buffer preparations were identical to the ones used for [68 Ga]Ga-PSMA (sodium ascorbate 0.4 M), with the addition of 300 µg ascorbic acid (0.3 mg in 0.1 mL water). 1.1 mL eluent was used for the extraction of the SCX cartridge. All synthesis reagents were contained in the reagent set except the peptide and ascorbic acid. From the reagent set, 50 mL of 0.9% NaCl was connected to the cassette on the designated spike. FAPI-46 was prepared in 10 µg, 20 µg, 30 µg, 40 µg and 50 µg aliquots from the stock solution (1 mg / mL). Meanwhile, the preparation of the acetate buffer solution was performed according to the user manual using EZ 102. The reaction mixture contains 0.4 mL of the final buffer solution, 0.1 mL of ascorbic acid and an aliquot of FAPI-46, which was loaded into the reaction vial. 3 mL of eluent was added to the eluent vial prior to synthesis.

Preparation for labeling of FAPI-46 with gallium-68 on ML eazy

First, the commercially available fully automated synthesis platform ML eazy was equipped with a disposable single-use cassette (PSMA cassette). FAPI-46 was prepared in $10-50\,\mu g$ aliquots from the stock solution (1 mg / mL). The buffer preparations were identical to the ones used for [68 Ga]Ga-PSMA, with the addition of 300 μg ascorbic acid (0.3 mg in 0.1 mL water). All synthesis reagents were contained in the reagent set except the peptide and ascorbic acid. From the reagent set, 6.0 mL of 0.9% NaCl was filled into the cassette in the designated vial with blue cap. After that, 1.1 mL of the eluent was filled into the designated vial with red cap. The reaction mixture contains 0.4 mL of the final buffer solution, 0.1 mL of ascorbic acid and an aliquot of FAPI-46, which was loaded into the reaction vial (green cap).

Labeling of FAPI-46 with 68GaCl₃

The synthesis was performed fully-automated without any user interaction. $^{68}\text{Ga}^{3^+}$ obtained from a 1.850 MBq $^{68}\text{Ge}/^{68}\text{Ga}$ generator (GalliaPharm*) with TiO_2 matrix, was eluted with 0.1 N HCl. The generator eluate was pre-concentrated on a strong cation exchange (SCX) cartridge. [$^{68}\text{GaCl}_4$] was recovered from the SCX cartridge by 1.1 mL eluent (5 M NaCl/HCl(0.1 M)) (Schultz et al. 2013). After the $^{68}\text{Ga-activity}$ was transferred to the reaction vial, it was heated to 95 °C. After 10 min at 95 °C (or 98 °C on ML

eazy), the reaction mixture was cooled down by adding 3.0 mL and 6.0 mL of 0.9% NaCl for MLPT and EAZY, respectively. The crude reaction solution was subsequently transferred through the CM cartridge and a sterile filter into the product vial. The reaction vessel was rinsed with 3.0 mL 0.9% NaCl in case of the MLPT synthesis. A sample for quality control was taken.

Quality control

After synthesis, the product vial was removed from the MLPT or ML eazy module and further evaluated for quality control determining the following parameters: total product activity, $^{68}\text{Ga}^{3+}$ -identity via half-life time, chemical purity (pH, sterility and endotoxins) and radiochemical purity (high-pressure-liquid chromatography). The stability of the product at room temperature was monitored by *radio*-HPLC for 3 h. *Radio*-HPLC was performed with the following method. A: acetonitrile + 0.1% TFA; B: water + 0.1% TFA, gradient: 0–8 min A: 0–50%; 8–10 min A: 50%; 10–12 min A: 50–0%. The limit of detection of the radio-HPLC is 10 kBq/ 20 μ L injected volume and the recovery is approximately 80%. For iTLC, ammonium acetate/methanol (1/1) and citrate buffer were used as mobile phases and iTLC-SG strips as stationary phase. The sterility tests were carried out as described in the European Pharmacopoeia by the Institut für Hygiene, Charité.

Abbreviations

FAP: Fibroblast activation protein; FAPI: Fibroblast activation protein inhibitor; n.d.c.: Non decay corrected; RCY: Radiochemical yield; MLPT: Modular Lab PharmTracer; ML eazy: Modular Lab eazy; PET: Positron Emission Tomography; ⁶⁸Ga^{3 +}: ⁶⁸Gallium³⁺ ion; ⁶⁸Ge^{4 +}: ⁶⁸Germanium⁴⁺ ion; HPLC: High-pressure-liquid-chromatography; TFA: Trifluoro acetic acid

Acknowledgements

We thank Eckert & Ziegler Eurotope GmbH for technical support. We thank SOFIE Biosciences for supplying the FAPI material and for technical support. We acknowledge support from the German Research Foundation (DFG) and the Open Access Publication Funds of Charité – Universitätsmedizin Berlin.

Authors' contributions

SS designed, performed and analyzed the experiments. MB and SP repeated the experiments in order to proof its repeatability. WB evaluated the results concerning clinical applications. All authors contributed to, read and approved the final manuscript.

Funding

FAPI-46 and the reference standard [natGa]Ga-FAPI-46 were made available free of costs from SOFIE Biosciences. Modular Lab eazy (ML eazy) was made available by Eckert & Ziegler within a research collaboration agreement with the Charité. Open Access funding enabled and organized by Projekt DEAL.

Availability of data and materials

All data generated or analyzed during this study are included in this published article.

Ethics approval and consent to participate

Not applicable.

Consent for publication

Not applicable.

Competing interests

The authors declare that they have no competing interests.

Received: 21 August 2020 Accepted: 20 November 2020 Published online: 17 December 2020

References

Boschi S, Lodi F, Malizia C, Cicoria G, Marengo M. Automation synthesis modules review. Appl Radiat Isot. 2013;76:38–45. Brennen WN, Isaacs JT, Denmeade SR. Rationale behind targeting fibroblast activation protein-expressing carcinoma-associated fibroblasts as a novel chemotherapeutic strategy. Mol Cancer Ther. 2012;11(2):257–66.

- Giesel F, Adeberg S, Syed M, Lindner T, Jimenez LD, Mavriopoulou E, Staudinger F, Tonndorf-Martini E, Regnery S, Rieken S, ElShafie R, Rohrich M, Flechsig P, Kluge A, Altmann A, Debus J, Haberkorn UA, Kratochwil C. FAPI-74 PET/CT using either (18)F-AIF or cold-kit (68)Ga-labeling: biodistribution, radiation Dosimetry and tumor delineation in lung Cancer patients. J Nucl Med. 2020;59(9):1415-22.
- Jansen K, Heirbaut L, Cheng JD, Joossens J, Ryabtsova O, Cos P, Maes L, Lambeir AM, De Meester I, Augustyns K, Van der Veken P. Selective inhibitors of fibroblast activation protein (FAP) with a (4-Quinolinoyl)-glycyl-2-cyanopyrrolidine scaffold. ACS Med Chem Lett. 2013;4(5):491–6.
- Lindner T, Loktev A, Altmann A, Giesel F, Kratochwil C, Debus J, Jager D, Mier W, Haberkorn U. Development of Quinoline-based Theranostic ligands for the targeting of fibroblast activation protein. J Nucl Med. 2018;59(9):1415–22.
- Loktev A, Lindner T, Burger EM, Altmann A, Giesel F, Kratochwil C, Debus J, Marme F, Jager D, Mier W, Haberkorn U. Development of fibroblast activation protein-targeted radiotracers with improved tumor retention. J Nucl Med. 2019; 60(10):1421–9.
- Meyer C, Dahlbom M, Lindner T, Vauclin S, Mona C, Slavik R, Czernin J, Haberkorn U, Calais J. Radiation dosimetry and biodistribution of (68)Ga-FAPI-46 PET imaging in cancer patients. J Nucl Med. 2019.
- Schultz MK, Mueller D, Baum RP, Leonard Watkins G, Breeman WA. A new automated NaCl based robust method for routine production of gallium-68 labeled peptides. Appl Radiat Isot. 2013;76:46–54.
- Velikyan I. 68Ga-based radiopharmaceuticals: production and application relationship. Molecules. 2015;20(7):12913-43.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Submit your manuscript to a SpringerOpen journal and benefit from:

- ► Convenient online submission
- ► Rigorous peer review
- ▶ Open access: articles freely available online
- ► High visibility within the field
- ► Retaining the copyright to your article

Submit your next manuscript at ▶ springeropen.com